

Test Definition: FPCEL

Elastase, Pancreatic, Serum

Overview

Method Name

Direct ELISA

NY State Available

No

Specimen

Specimen Type

Serum

Specimen Required

Patient Preparation:

- 1. Fasting: 10 to 12 hours, recommended
- 2. Medications that affect pancreatic activity should be discontinued, if possible, for at least 48 hours prior to specimen collection.

Collection Container/Tube:

Preferred: Red top **Acceptable**: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL serum

Collection Instructions:

- 1. As soon as possible after collection, centrifuge and aliquot 3 mL of serum into a plastic vial.
- 2. Freeze immediately.
- 3. Send frozen.

Specimen Minimum Volume

Serum: 1 mL

Reject Due To

Gross	Reject
hemolysis	
Gross lipemia	Reject
Gross icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
• • • • • • • • • • • • • • • • • • • •	-		-



Test Definition: FPCEL

Elastase, Pancreatic, Serum

Serum	Frozen (preferred)	365 days	
	Ambient	72 hours	
	Refrigerated	7 days	

Clinical & Interpretive

Clinical Information

Refer to www.interscienceinstitute.com/individual-assays/

Reference Values

Normal pancreatic exocrine function: Less than 3.5 ng/mL

No pediatric reference ranges are available for this test.

Performance

PDF Report

Referral

Day(s) Performed

Monday through Friday

Report Available

12 to 14 days

Performing Laboratory Location

Inter Science Institute

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test has not been cleared or approved by the US Food and Drug Administration.

This test was developed and its performance characteristics determined by Inter Science Institute. Values obtained with different methods, laboratories, or kits cannot be used interchangeably with the results on this report. The results cannot be interpreted as absolute evidence of the presence or absence of malignant disease.



Test Definition: FPCEL

Elastase, Pancreatic, Serum

CPT Code Information

83520

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPCEL	Pancreatic Elastase-1	Not Provided
Result ID	Test Result Name	Result LOINC® Value