

Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC) Confirmation, Random, Urine

## Overview

## **Useful For**

Detection and confirmation of drug use of cannabis/marijuana and to specifically identify and quantify delta-8-carboxy tetrahydrocannabinol (THC-COOH) and delta-9-THC-COOH

## **Special Instructions**

• Clinical Toxicology CPT Code Client Guidance

### **Method Name**

Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS)

## **NY State Available**

Yes

# Specimen

# **Specimen Type**

Urine

## **Ordering Guidance**

For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order THCX / Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC) Confirmation, Chain of Custody, Random, Urine

Additional drug panels and specific requests are available. Call 800-533-1710 or 507-266-5700.

# **Additional Testing Requirements**

If urine creatinine is required or adulteration of the sample is suspected, order, ADULT / Adulterants Survey, Random, Urine.

# **Specimen Required**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container Tube: Plastic urine container

Submission Container/Tube: Plastic vial

**Specimen Volume:** 3 mL Collection Instructions:

- 1. Collect a random urine specimen.
- 2. No preservative.

**Additional Information:** 



Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC) Confirmation, Random, Urine

- 1. No specimen substitutions.
- STAT requests are not accepted for this test.

#### **Forms**

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

- -General Request (T239)
- -Therapeutics Test Request (T831)

## Specimen Minimum Volume

0.5 mL

## Reject Due To

Gross	ОК
hemolysis	
Gross icterus	Reject

# **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	14 days	

# **Clinical & Interpretive**

#### **Clinical Information**

There are over 100 different cannabinoids in cannabis/marijuana. The main psychoactive cannabinoid is delta-9-tetrahydrocannabinol (delta-9-THC), which is the active agent of the popularly abused street drug, cannabis/marijuana. Delta-8 tetrahydrocannabinol (delta-8 THC) is another psychoactive substance found in the *Cannabis sativa* plant, of which cannabis/marijuana and hemp are 2 varieties. Delta-8 THC is one of over 100 cannabinoids produced naturally by the cannabis plant but is not typically found in significant amounts in the plant itself. As a result, concentrated amounts of delta-8 THC are typically manufactured from hemp-derived cannabidiol.

Following consumption of cannabis/marijuana, delta-9-THC metabolizes to a variety of inactive products, one of them being the carboxy metabolite (delta-9-THC-COOH). In almost all medico-legal cases or when the patient adamantly denies cannabis/marijuana use and the immunoassay test is positive, confirmation of the result by a definitive test is required. This test is a definitive, confirmatory test using liquid chromatography tandem mass spectrometry to identify and quantify delta-8-THC-COOH and delta-9 THC-COOH.

## **Reference Values**

Not Detected (Positive results are reported with a quantitative result.)



Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC) Confirmation, Random, Urine

Cutoff concentration by liquid chromatography tandem mass spectrometry:

Delta-8-Carboxy-Tetrahydrocannabinol (THC): 5.0 ng/mL

Delta-9-Carboxy-Tetrahydrocannabinol (THC): 5.0 ng/mL

# Interpretation

The presence of delta-8 and/or delta-9 carboxy tetrahydrocannabinol (THC-COOH) in urine is a strong indicator that the patient has used cannabis/marijuana.

THC-COOH has a long half-life and can be detected in urine for more than 7 days after a single use.

Chronic use causes accumulation of THC and THC-COOH in adipose tissue, such that it is excreted into the urine for as long as 30 to 60 days from the time chronic use is halted.

### **Cautions**

No significant cautionary statements

### Clinical Reference

- 1. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 12th ed. Biomedical Publications; 2020
- 2. Langman LJ, Bechtel LK, Holstege CP. Clinical toxicology. In: Rifai N, Chiu RWK, Young I, Burnham CAD, Wittwer CT, eds. Tietz Textbook of Laboratory Medicine. 7th ed. Elsevier; 2023:chap 43

## **Performance**

## **Method Description**

Carboxy tetrahydrocannabinol (THC-COOH) is extracted from urine by making the urine alkaline to hydrolyze THC-COOH glucuronide. The hydrolyzed samples are diluted for analysis by liquid chromatography tandem mass spectrometry using selected ion monitoring. (Unpublished Mayo method)

# **PDF Report**

No

# Day(s) Performed

Monday through Sunday

### Report Available

3 to 5 days

# **Specimen Retention Time**

2 weeks

## **Performing Laboratory Location**



Delta-8 and Delta-9-Carboxy-Tetrahydrocannabinol (THC)
Confirmation, Random, Urine

Mayo Clinic Laboratories - Rochester Superior Drive

# **Fees & Codes**

# **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

## **Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

## **CPT Code Information**

G0480

80349 (if appropriate for select payers)

Clinical Toxicology CPT Code Client Guidance

# **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
THCU	Carboxy-THC Confirmation, U	20521-1

Result ID	Test Result Name	Result LOINC® Value
2497	Delta-9	20521-1
	Carboxy-Tetrahydrocannabinol by	
	LC-MS/MS	
21186	Carboxy-THC Interpretation	69050-3
618770	Delta-8	20521-1
	Carboxy-Tetrahydrocannabinol by	
	LC-MS/MS	