

Test Definition: FLIPR

Lipase, Random Urine

Overview

Method Name

Spectrophotometry

NY State Available

Yes

Specimen

Specimen Type

Urine

Specimen Required

Container/Tube: Plastic urine container

Specimen Volume: 2 mL Collection Instructions:

- 1. Collect random urine without preservative.
- 2. Ship 2 mL urine in a sterile screw capped plastic container ship ambient.

Specimen Minimum Volume

1 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	Preservatives

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Ambient (preferred)	5 days	
	Refrigerated	5 days	
	Frozen	21 days	

Clinical & Interpretive



Test Definition: FLIPR

Lipase, Random Urine

Reference Values

Adult: < or = 4 U/L

Performance

PDF Report

No

Day(s) Performed

Tuesday - Friday, Sunday

Report Available

1 to 7 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

CPT Code Information

83690

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FLIPR	Lipase, Urine	10888-6

Result ID	Test Result Name	Result LOINC® Value
Z1559	Lipase, Urine	10888-6