

---

## Human Leukocyte Antigens (HLA)-DR-DQ Disease Association Typing Low Resolution, Blood

**Test ID:** 2DIS

**Useful for:**

Determining class II human leukocyte antigens (HLA) antigens to identify potential disease association and as markers for drug hypersensitivity

**Method:**

Polymerase Chain Reaction (PCR)/Sequence-Specific Oligonucleotide Probes (SSO)

**Reference Values:**

Not applicable

**Specimen Requirements:**

**Container/Tube:** Yellow top (ACD Solution B)

**Specimen Volume:** 6 mL

**Collection Instructions:** Send in original tube. Do not transfer blood to other containers.

**Additional Information:** Specimen acceptability is based on extracted DNA concentration and not sample age.

**Minimum Volume:** 3 mL

**Specimen Stability Information:**

Specimen Type	Temperature	Time
Whole blood ACD	Ambient	
	Refrigerated	

**Cautions:**

No significant cautionary statements.

**CPT Code:**

81376 x4

**Day(s) Setup:** Monday through Friday 7:30 a.m.-5p.m.      **Analytic Time:** 5 days

**Questions**

Contact Richard Einerson, Laboratory Technologist Resource Coordinator at 800-533-1710.