

## Hemoglobin Electrophoresis Evaluation, Blood

**Test ID:** HBELP

**Note:** This test will replace HBELX / Hemoglobin Electrophoresis Cascade, Blood.

The following Reflex Tests will have a Reporting Name change:

- IEF – Isoelectric Focusing, B
- HPFH – Hb F Distribution, B
- SDEX – Sickie Solubility, B
- UNHB – Hb Stability, B

### Useful for:

Diagnosis and classification of hemoglobin disorders, including thalassemias and hemoglobin variants

### Profile Information:

Test ID	Reporting Name	Available Separately	Always Performed
HBELI	Hb Electrophoresis Interpretation	No	Yes
HGBCE	Hb Variant, A2 and F Quantitation,B	Yes	Yes
HPLC	HPLC Hb Variant, B	No	Yes

### Reflex Tests:

Test ID	Reporting Name	Available Separately	Always Performed
HPFH	Hb F Distribution, B	No	No
MASS	Hb Variant by Mass Spec, B	No	No
SDEX	Sickle Solubility, B	Yes	No
IEF	Isoelectric Focusing, B	No	No
UNHB	Hb Stability, B	No	No
ATHAL	Alpha-Globin Gene Analysis	Yes	No
WASQR	Alpha Globin Gene Sequencing, B	Yes (Order WASEQ)	No
WBSQR	Beta Globin Gene Sequencing, B	Yes (Order WBSEQ)	No
WBDDR	Beta Globin Cluster Locus Del/Dup,B	Yes (Order WBDD)	No
WGSQR	Gamma Globin Full Gene Sequencing	No	No
HBEL0	Hb Electrophoresis Summary Interp	No	No

## Methods:

Capillary Electrophoresis (HGBCE), Cation Exchange/High-Performance Liquid Chromatography (HPLC), Isoelectric Focusing (IEF), Mass Spectrometry (MASS), Flow Cytometry (HPFH), Isopropanol and Heat Stability (UNHB).

## Reference Values:

### HEMOGLOBIN A:

1-30 days: 5.9-77.2%  
1-2 months: 7.9-92.4%  
3-5 months: 54.7-97.1%  
6-8 months: 80.0-98.0%  
9-12 months: 86.2-98.0%  
13-17 months: 88.8-98.0%  
18-23 months: 90.4-98.0%  
> or =24 months: 95.8-98.0%

### HEMOGLOBIN A2:

1-30 days: 0.0-2.1%  
1-2 months: 0.0-2.6%  
3-5 months: 1.3-3.1%  
> or =6 months: 2.0-3.3%

### HEMOGLOBIN F:

1-30 days: 22.8-92.0%  
1-2 months: 7.6-89.8%  
3-5 months: 1.6-42.2%  
6-8 months: 0.0-16.7%  
9-12 months: 0.0-10.5%  
13-17 months: 0.0-7.9%  
18-23 months: 0.0-6.3%  
> or =24 months: 0.0-0.9%

### VARIANT 1

0.0

### VARIANT 2

0.0

### VARIANT 3

0.0

## Specimen Requirements:

### Container/Tube:

**Preferred:** Lavender top (EDTA)

**Acceptable:** ACD (solution B), green top (sodium heparin)

**Specimen Volume:** 10 mL

**Collection Instructions:** Send specimen in original tube. **Do not aliquot.**

**Minimum Volume:** 1 mL (this volume will limit reflex testing possibilities)  
3 mL if multiplex ligation-dependent probe amplification is desired

## Specimen Stability Information:

Specimen Type	Temperature	Time
Whole Blood EDTA	Refrigerated	7 days

## Cautions:

Some hemoglobin disorders and variants are not detected by the screening methods including common alpha thalassemia conditions and require further reflex testing to identify. If a family history of a known hemoglobin disorder, prior therapy for a hemoglobin disorder, or otherwise unexplained lifelong/familial symptoms such as hemolysis, microcytosis, erythrocytosis/polycythemia, cyanosis, or hypoxia are present, this should be clearly communicated to the laboratory so appropriate reflex testing can be added, see [Metabolic Hematology Patient Information](#) (T810).

Recent transfusion may mask protein results including hemoglobin electrophoresis, hereditary persistence of hemoglobin F (HPFH) by flow cytometry, stability studies, and sickle solubility studies depending on percentage of transfused cells present.

Some hemoglobin variants can originate from the donor blood product and not from the tested recipient. These are typically found in low percentage.

If the patient has undergone a bone marrow transplant, the results may show atypical results and should be interpreted in the context of clinical information.

Some therapies cause artefactual effects in protein studies, including hydroxyurea and decitabine (increased Hb F levels), Voxelotor (artefactual peaks) and gene therapy (alternate protein detection, Beta T87Q, by mass spectrometry). Clear communication of prior therapy is strongly recommended.

**CPT Code:**

83020-Quantitation by electrophoresis

83021-Quantitation by HPLC

82664-Electrophoresis, not elsewhere specified (if appropriate)

83068 (if appropriate)

83789 (if appropriate)

88184 (if appropriate)

**Day(s) Setup:** Monday through Saturday

**Analytic Time:** 2 days

**Questions**

Contact Amy Nelson, Laboratory Technologist Resource Coordinator at 800-533-1710.