

## Fentanyl with Metabolite Confirmation, Random, Urine

**Test ID:** FENTU

**Explanation:** A method enhancement will result in the following test changes.

Current Specimen Required
<b>Supplies:</b> Urine Tubes, 10 mL (T068) <b>Collection Container/Tube:</b> Clean, plastic urine collection container <b>Submission Container/Tube:</b> Plastic, 10-mL urine tube <b>Specimen Volume:</b> 10 mL <b>Specimen Minimum Volume:</b> 2.1 mL <b>Collection Instructions:</b> 1. Collect a random urine specimen. 2. No preservative. <b>Additional Information:</b> 1. No specimen substitutions. 2. No STATS are accepted for this procedure. 3. Submitting less than 10 mL will compromise our ability to perform all necessary testing.

New Specimen Required
<b>Supplies:</b> Sarstedt Aliquot Tube, 5 mL (T914) <b>Collection Container/Tube:</b> Clean, plastic urine collection container <b>Submission Container/Tube:</b> Plastic, 5-mL tube <b>Specimen Volume:</b> 5 mL <b>Specimen Minimum Volume:</b> 1 mL <b>Collection Instructions:</b> 1. Collect a random urine specimen. 2. No preservative. <b>Additional Information:</b> 1. No specimen substitutions. 2. No STATS are accepted for this procedure. 3. Submitting less than 5 mL may compromise our ability to perform all necessary testing.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	14 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	21 days
	Frozen	28 days

Current Result Codes	
Result ID	Reporting Name
45530	Fentanyl
45531	Norfentanyl
45532	Fentanyl Interpretation

New Result Codes	
Result ID	Reporting Name
45530	Fentanyl by LC-MS/MS
45531	Norfentanyl by LC-MS/MS
45532	Fentanyl Interpretation

### Questions

Contact Brandon DeBoom, Laboratory Technologist Resource Coordinator at 800-533-1710.