

Angiotensin I, Plasma

Explanation: The following new orderable, referred to Inter Science Institute, will be available on 11/12/2024.

Test ID: FANG1

Methods:

Direct ELISA

Reference Values:

Up to 25 pg/mL

Specimen Requirements:

Patient preparation: Patient should be on a normal sodium diet, 110 mEq of sodium. Patient should be in a recumbent posture for at least 30 minutes prior to collection of specimen. Diuretics, mineralocorticoids, glucocorticoids, estrogens, oral contraceptives, and ACTH medications and sodium, potassium, and posture all affect Angiotensin levels.

Specimen Type: Plasma

Collection Container/Tube: EDTA

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood in a lavender-top (EDTA) tube(s)
2. Centrifuge and aliquot 3 mL of plasma into a plastic vial.
3. Send frozen.

Minimum Volume: 1 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Plasma EDTA	Frozen	30 days

Fee:
\$180.00

CPT Code:
83520

Day(s) Setup: Monday through Friday

Report Available: 12 to 14 days

Questions

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.