

## Titanium/Creatinine Ratio, Random, Urine

**Test ID:** TIUCR

**Explanation:**

TIUCR will become obsolete on the effective date.

**Recommended Alternative Test:**

## Titanium, Urine

**Test ID:** ZW86

**Referral Lab Code:** 4486U

**Methodology:**

Inductively Coupled Plasma/Mass Spectrometry (ICP/MS)

**Reference Values:**

Interpretation provided on report.

**Specimen Requirements:**

**Container/Tube:** Plastic container

**Preferred:** Plastic container (Acid washed or Trace metal-free)

**Acceptable:** Plastic container (Acid washed or Trace metal-free)

**Specimen Volume:** 2 mL

**Collection Instructions:** Collect urine directly into a labeled trace-metal free (acid-washed, deionized water-rinsed, air-dried), polyethylene bottle with a plastic lined screw cap or an acid washed collection container.

PLEASE NOTE: Avoid exposure to gadolinium-based contrast media for 48 hours prior to sample collection

**Minimum Volume:** 0.6 mL

**Specimen Stability Information:**

<b>Specimen Type</b>	<b>Temperature</b>	<b>Time</b>
Urine	Frozen	30 days

**CPT Code:**83018

**Miscellaneous Referral Test Fee:** \$212.00

**Report Available:** 13 days

**Questions**

Contact Rebekah Walsh, Laboratory Resource Coordinator at 800-533-1710.