

Overview

Method Name
Enzyme-Linked Immunosorbent Assay (ELISA)

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
Serum: 0.4 mL
Note: This volume does not allow for repeat testing.

Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Specimen other than serum	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	

	Ambient	14 days	
	Frozen	14 days	

Clinical & Interpretive

Clinical Information

Refer to www.labcorp.com/test-menu/search

Reference Values

0.0-5.0 units/mL

Cautions

This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H or coenzyme R). It is recommended to ask all patients who may be indicated for this test about biotin supplementation. Patients should be cautioned to stop biotin consumption at least 72 hours prior to the collection of a sample.

Performance

PDF Report

No

Day(s) Performed

Monday through Saturday

Report Available

5 to 8 days

Performing Laboratory Location

LabCorp Burlington

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86341

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FGADA	GAD-65 Autoantibody	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FGADA	GAD-65 Autoantibody	56540-8