MAYO CLINIC LABORATORIES

Test Definition: FDA1S

Drugs of Abuse (10 Panel) and Alcohol Screen, Serum

Overview

Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FAMCS	Amphetamines Confirmation, S	No	No
FBCS	Benzodiazepines Confirmation, S	No	No
FCANS	Cannabinoids Confirmation, S	No	No
FCMC	Cocaine and Metabolites Confirm, S	No	No
FETCS	Ethanol Confirmation, S	No	No
FMMCS	Methadone and Metabolite Confirm, S	No	No
FOFU	Opiates-Free Unconjugated Conf, S	No	No
FPCCS	Phencyclidine Confirmation, S	No	No
FBCFS	Barbiturates Confirmation, S	No	No

Testing Algorithm

If appropriate, confirmation testing will be performed at an additional charge.

Method Name

Enzyme-Linked Immunosorbent Assay (ELISA)
Headspace Gas Chromatography (GC)

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Test Definition: FDA1S

Drugs of Abuse (10 Panel) and Alcohol Screen, Serum

*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 4 mL Serum
Collection Instructions:
1. Collect specimen using alcohol-free skin preparation.
2. Promptly centrifuge and aliquot serum into a plastic vial.
3. Send refrigerate.

Specimen Minimum Volume
Serum: 2 mL

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	14 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog.

Reference Values
Reporting limit determined each analysis.

Ethanol: None detected
Opiates: None detected
Cocaine/Metabolites: None detected
Benzodiazepines: None detected
Cannabinoids: None detected
Amphetamines: None detected
Barbiturates: None detected
Methadone/Metabolite: None detected
Phencyclidine: None detected
Methamphetamine/MDMA: None detected
Oxycodone/Oxymorphone: None detected

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
6 to 20 days

Performing Laboratory Location
NMS Labs

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information
80307x2

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDA1S	Drug of Abuse 10 and Alcohol Scrn,S	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5535	Ethanol	20470-1
Z5536	Opiates	8219-8
Z5537	Cocaine/Metabolites	73823-7
Z5538	Benzodiazepines	46976-7
Z5539	Cannabinoids	8172-9
Z5540	Amphetamines	8149-7

Z5541	Barbiturates	20421-4
Z5542	Methadone/Metabolite	59705-4
Z5543	Phencyclidine	8236-2
Z5544	Methamphetamine/MDMA	91034-9
Z5545	Oxycodone/Oxymorphone	77732-6