MAYO CLINIC LABORATORIES

Test Definition: FRT3S

T3 Reverse, LC/MS/MS

Overview

Method Name
Chromatography/Mass Spectrometry

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube: Red top or serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Allow blood to clot at room temperature for 30 minutes.
2. Within 24 hours of collection, centrifuge and aliquot serum immediately into a plastic vial.
3. Send frozen.

Specimen Minimum Volume
Serum: 0.3 mL

Reject Due To

Specimen received in serum gel separator (SST) tube	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	30 days	
	Ambient	7 days	
	Refrigerated	7 days	

Test Definition: FRT3S

T3 Reverse, LC/MS/MS

Clinical & Interpretive

Clinical Information

Refer to https://testdirectory.questdiagnostics.com/test/home

Reference Values

8-25 ng/dL

Performance

PDF Report

No

Day(s) Performed

Daily

Report Available

4 to 8 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed, and its analytical performance characteristics have been determined by Quest Diagnostics. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.

CPT Code Information

84482

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FRT3S	T3 Reverse	3052-8

MAYO CLINIC LABORATORIES

Test Definition: FRT3S

T3 Reverse, LC/MS/MS

Result ID	Test Result Name	Result LOINC® Value
FRT3S	T3 Reverse	3052-8