MAYO CLINIC LABORATORIES

Test Definition: FFLUR

Flurazepam (Dalmane) and Desalkylflurazepam

Overview

Method Name
Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
Submit only 1 of the following specimens:

Plasma
Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 2 mL sodium heparin plasma frozen in plastic vial.

Serum
Draw blood in a plain, red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 2 mL of serum frozen in plastic vial.

Specimen Minimum Volume
0.3 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Frozen	180 days	

Clinical & Interpretive

Reference Values

Flurazepam:
Reference Range: 0 - 30 ng/mL

Desalkylflurazepam:
Reference Range: 30 - 150 ng/mL

Flurazepam + Desalkyflurazepam:
Reference Range: 30 - 180 ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 9 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

80346

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFLUR	Flurazepam (Dalmane)	73826-0

MAYO CLINIC LABORATORIES

Test Definition: FFLUR

Flurazepam (Dalmane) and Desalkylflurazepam

Result ID	Test Result Name	Result LOINC® Value
Z1133	Flurazepam	46979-1
Z1160	Desalkylflurazepam	46980-9
Z1161	Flurazepam + Desalkylflurazepam	73826-0