



PATIENT NAME TESTRNV, IMPLEMENTATION				ORDER NUMBER M522000364
PATIENT ID X100404241	DATE OF BIRTH 02/19/1968	AGE 54 Y	SEX Male	REQUESTED BY ATLAS TEST
COLLECTED 12/22/2022, 6:33 AM	RECEIVED 12/22/2022, 2:02 PM	REPORTED 1/3/2023, 10:33 AM		
The collected, received, and reported dates and times on the report are in the time zone of the performing location. 7028846 MCL RochesterCampus Rochester MN 55901				CLIENT ORDER NUMBER X100404241 CLIENT MRN 321

TEST DESCRIPTION

Evaluation of the ATP7B gene associated with Wilson disease

SPECIMEN

WB Whole Blood

RESULT SUMMARY

Pathogenic Variants Detected

RESULT

Gene (Transcript)	Variant	Zygosity	Classification
ATP7B (NM_000053.4)	c.3207C>A p.His1069Gln (p.H1069Q) Chr13(GRCh37):g.52518281G>T	homozygous	PATHOGENIC

The following homozygous PATHOGENIC variant was detected:

ATP7B (NM_000053.4), Chr13(GRCh37):g.52518281G>T, c.3207C>A, p.His1069Gln (p.H1069Q)

INTERPRETATION

ATP7B c.3207C>A (p.His1069Gln), PATHOGENIC

The homozygous c.3207C>A (p.His1069Gln) missense variant in the ATP7B gene (MIM:606882) is an established pathogenic variant. Pathogenic variants in the ATP7B gene have been associated with autosomal recessive Wilson disease (1). This result is supportive of a diagnosis of Wilson disease for this individual. Clinical correlation is recommended.

This result should be interpreted in the context of clinical findings, family history, and other laboratory testing.

Consultation with a genetics professional is recommended for interpretation of this result and to determine whether reproductive risk assessment and familial testing may be of benefit to this family. Genetic testing for family members is available by ordering FMTT / Familial Mutation, Targeted Testing for the specific variant(s) detected. Please contact the laboratory at 1-800-533-1710 or the online test catalog at www.mayocliniclabs.com for information about FMTT.

REFERENCES:

1) GeneReviews: Wilson Disease (<https://www.ncbi.nlm.nih.gov/books/NBK1512/>) (PMID: 20301685)

METHOD

Next generation sequencing (NGS) and/or Sanger sequencing was performed to test for the presence of variants in coding regions and intron/exon boundaries of the gene analyzed. The human genome reference GRCh37/hg19 build was used for sequence read alignment. At least 99% of the bases are covered at a read depth >30X. Sensitivity is estimated at >99% for single nucleotide variants, >94% for indels up to 39 base pairs, >95% for deletions up to 75 base pairs and insertions up to 47 base pairs. NGS and/or a PCR-based quantitative method was performed to test for the presence of

PATIENT NAME	TESTRNV, IMPLEMENTATION	22-49998
--------------	-------------------------	----------



1-800-533-1710

ATP7B Full Gene Analysis

WNDZ

PATIENT NAME TESTRNV, IMPLEMENTATION				ORDER NUMBER M522000364
PATIENT ID X100404241	DATE OF BIRTH 02/19/1968	AGE 54 Y	SEX Male	REQUESTED BY ATLAS TEST
COLLECTED 12/22/2022, 6:33 AM	RECEIVED 12/22/2022, 2:02 PM	REPORTED 1/3/2023, 10:33 AM		
The collected, received, and reported dates and times on the report are in the time zone of the performing location. 7028846 MCL RochesterCampus Rochester MN 55901				CLIENT ORDER NUMBER X100404241 CLIENT MRN 321

deletions and duplications in the gene analyzed. See the Genes Analyzed field for a list of gene(s) tested.

There may be regions of genes that cannot be effectively evaluated for sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high GC content, and repetitive sequences. Confirmation of select reportable variants was performed by alternate methodologies based on internal laboratory criteria. See www.mayocliniclabs.com (TEST ID WNDZ) for details regarding genes with regions not routinely covered.

GENES ANALYZED

ATP7B

DISCLAIMER**Clinical Correlations**

An online research opportunity called GenomeConnect (genomeconnect.org), a project of ClinGen, is available for the recipient of this genetic test. This patient registry collects de-identified genetic and health information to advance the knowledge of genetic variants. Mayo Clinic is a collaborator of ClinGen. This may not be applicable for all tests.

If testing was performed because of a clinically significant family history it is often useful to first test an affected family member. Detection of a reportable variant(s) in an affected family member would allow for more informative testing of at risk individuals.

To discuss the availability of further testing options or for assistance in the interpretation of these results, Mayo Clinic Laboratory genetic counselors can be contacted at 1-800-533-1710.

Technical Limitations

Next generation sequencing may not detect all types of genomic variants. In rare cases, false negative or false positive results may occur. The depth of coverage may be variable for some target regions, but assay performance below the minimum acceptable criteria or for failed regions will be noted. Given these limitations, negative results do not rule out the diagnosis of a genetic disorder. If a specific clinical disorder is suspected, evaluation by alternative methods can be considered.

There may be regions of genes that cannot be effectively evaluated for sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high GC content, and repetitive sequences. Confirmation of select reportable variants was performed by alternate methodologies based on internal laboratory criteria.

Additionally, low level mosaic variants may not be detected.

This test is not designed to differentiate between somatic and germline variants. If there is a possibility that any detected variant is somatic, additional testing may be necessary to clarify the significance of results.

Reclassification of Variants Policy

See www.mayocliniclabs.com (TEST ID WNDZ) for information regarding the laboratory's policy for reclassification of variants.



PATIENT NAME TESTRNV, IMPLEMENTATION					ORDER NUMBER M522000364
PATIENT ID X100404241	DATE OF BIRTH 02/19/1968	AGE 54 Y	SEX Male	REQUESTED BY ATLAS TEST	
COLLECTED 12/22/2022, 6:33 AM	RECEIVED 12/22/2022, 2:02 PM	REPORTED 1/3/2023, 10:33 AM			
The collected, received, and reported dates and times on the report are in the time zone of the performing location.				CLIENT ORDER NUMBER X100404241	
7028846 MCL RochesterCampus Rochester MN 55901				CLIENT MRN 321	

Variant Evaluation

Variant curation is performed using published ACMG-AMP recommendations as a guideline. Other gene-specific guidelines may also be considered. Variants classified as benign or likely benign are not reported.

Results from in silico evaluation tools may change over time and should be interpreted with caution and professional clinical judgment.

TEST CLASSIFICATION

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

RELEASED BY

Devin Oglesbee, Ph.D.

Code : MCR Laboratory : Mayo Clinic Laboratories - Rochester Main Campus
Lab Director : WILLIAM G MORICE, II MD, PhD CLIA Certificate : 24D0404292

Address : 200 FIRST STREET SW
ROCHESTER MN 55905