

## Overview

### Method Name

Enzyme Immunoassay

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

### Reject Due To

Hemolysis	NA
Lipemia	Mild OK; Gross reject
Icterus	NA
Other	NA

### Specimen Minimum Volume

0.5 mL

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

## Clinical & Interpretive

### Reference Values

Negative

The Phadia Phadiatop test is an allergy screening test with excellent sensitivity and specificity for inhalant allergy. It uses an ImmunoCAP with a balanced mixture of representative allergens, including grasses, trees, weeds, cat, dog, mites and molds. A positive result indicates that the patient is allergic to one or more of these allergens; a negative indicates the patient is not allergic to inhalant allergens. Note that the test does not assess a patient's sensitivity to food, drug, chemical or certain unusual or rare allergens.

## Performance

### PDF Report

No

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Test Classification

### CPT Code Information

86005