

**Overview****Method Name**

Enzyme Immunoassay

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Hemolysis	NA
Lipemia	Mild OK; Gross reject
Icterus	NA
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Clinical & Interpretive****Reference Values**

Negative

The Phadia Phadiatop test is an allergy screening test with excellent sensitivity and specificity for inhalant allergy. It uses an ImmunoCAP with a balanced mixture of representative allergens, including grasses, trees, weeds, cat, dog, mites and molds. A positive result indicates that the patient is allergic to one or more of these allergens; a negative indicates the patient is not allergic to inhalant allergens. Note that the test does not assess a patient's sensitivity to food, drug, chemical or certain unusual or rare allergens.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Friday

### Report Available

1 to 6 days

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### CPT Code Information

86005

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPHAS	Phadiatop (Allergy Screen)	57036-6

Result ID	Test Result Name	Result LOINC® Value
FPHAS	Phadiatop (Allergy Screen)	57036-6