

Overview**Method Name**

Gel Diffusion

NY State Available

Yes

Specimen**Specimen Type**

Serum

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Clinical & Interpretive**Reference Values**

Canary Droppings

Negative

Chicken Serum	Negative
Cockatiel Droppings	Negative
Finch Droppings	Negative
Parakeet Droppings	Negative
Parakeet Serum	Negative
Parrot Droppings	Negative
Parrot Serum	Negative
Pigeon/Dove Droppings	Negative
Pigeon/Dove Serum	Negative

The gel diffusion method was used to test this patient's serum for the presence of precipitating antibodies (IgG) to the antigens indicated. These antibodies are serological markers for exposure and immunological sensitization. The clinical significance varies, depending on the history and symptoms.

Performance

PDF Report

No

Day(s) Performed

Tuesday, Friday

Report Available

3 to 11 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86331 x 10

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBFPI	Bird Fancier's Precipitin Panel I	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z3892	Parakeet Droppings Gel Diffusion	24503-5
Z3894	Parrott Droppings Gel Diffusion	24501-9
Z3890	Cockatiel Droppings Gel Diffusion	23843-6
Z3891	Finch Droppings Gel Diffusion	Not Provided
Z3893	Parakeet Serum Gel Diffusion	24502-7
Z3888	Canary Droppings Gel Diffusion	Not Provided
Z3889	Chicken Serum Gel Diffusion	24498-8
Z3895	Parrott Serum Gel Diffusion	24500-1
Z3897	Pigeon/Dove Serum Gel Diffusion	6733-0
Z3896	Pigeon/Dove Droppings Gel Diffusion	24505-0