

Overview

Method Name

Hemolysis

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Specimen Type: Serum

Container/Tube: Red

Specimen Volume: 3 mL

Draw blood in a plain red-top tube. Maintain specimen at 37 degrees C until serum is separated from cells. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	Separator or gel tube

Specimen Minimum Volume

2 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Ambient		

Clinical & Interpretive**Reference Values**

Negative

Performance**PDF Report**

No

Performing Laboratory Location

ARUP Laboratories

Fees & Codes**Test Classification****CPT Code Information**

86940, 86941

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FDLS	Donath Landsteiner	18287-3

Result ID	Reporting Name	LOINC®
FDLS	Donath Landsteiner	18287-3