

Overview
Testing Algorithm

If the GHB Screen, UR is non-negative, then the Gamma-Hydroxybutyric Acid, CF, UR (FGHAC) will be performed at an additional charge.

Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FGHAC	Gamma-Hydroxybutyric Acid, CF, UR	No	No

Method Name

Liquid Chromatography with Tandem Mass Spectrometry (LC/MS/MS)

Gas Chromatography/Mass Spectrometry (GC/MS) (if appropriate)

NY State Available

Yes

Specimen
Specimen Type

Urine

Specimen Required

Collect 10 mL random urine without preservatives. Ship refrigerated in a plastic container.

Reject Due To

Hemolysis NA

Lipemia NA

Icterus NA

Other NA

Specimen Minimum Volume

1.2 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient		

Clinical & Interpretive

Reference Values

Reference Range: Negative

Screening threshold: 5.0 ug/mL

Performance**PDF Report**

No

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes**Test Classification****CPT Code Information**

80307