

**Overview**
**Testing Algorithm**

If the Cocaine Analysis, whole blood is positive, then the Cocaine, Confirmation (FCOKR) will be performed at an additional charge.

**Reflex Tests**

Test Id	Reporting Name	Available Separately	Always Performed
FCOKR	Cocaine, MS, WB/SP Rfx	No	No

**Method Name**

Immunoassay (IA)

Gas Chromatography with mass spectrometry (GC/MS); confirmation

**NY State Available**

Yes

**Specimen**
**Specimen Type**

Whole blood

**Specimen Required**

Draw blood in a purple-top (EDTA), green-top (sodium-heparin) or grey-top (NaF/oxalate) tube(s). Send 7 mL of whole blood refrigerated in a plastic vial.

**Reject Due To**

Hemolysis NA

Lipemia NA

Icterus NA

Other NA

**Specimen Minimum Volume**

3 mL

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Whole blood	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient		

**Clinical & Interpretive****Reference Values**

Cocaine Screen:

Cocaine and Metabolite, UA â€œ Negative; Cutoff: 25 ng/mL

Cocaine Confirmation:

Cocaine

Benzoylcegonine

Confirmation threshold: 10 ng/mL

**Performance****PDF Report**

No

**Performing Laboratory Location**

Medtox Laboratories, Inc.

**Fees & Codes****Test Classification**

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

**CPT Code Information**

80307

80353-if applicable

**LOINC® Information**

Test ID	Test Order Name	Order LOINC Value
FCOKE	Cocaine, Scr w/Conf, WB	16633-0

Result ID	Test Result Name	Result LOINC Value
Z4826	Cocaine and Metabolite	73823-7