

Overview

Method Name

Next Generation Sequencing

NY State Available

Yes

Specimen

Specimen Type

Whole Blood EDTA

Specimen Required

Draw blood in a lavender-top (EDTA) tube(s) and send 8 mL of whole blood. Ship Ambient.

Note: Collection date and informed consent are required.

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	Frozen specimens

Specimen Minimum Volume

6 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Ambient (preferred)	10 days	
	Refrigerated	10 days	

Clinical & Interpretive

Reference Values

A final report will be attached in MayoAccess.

Performance**PDF Report**

Referral

Performing Laboratory Location

Athena Diagnostics

Fees & Codes**Test Classification****CPT Code Information**

81406

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FSCN4	SCN4A (Myotonia) DNA Sequencing	Not Provided

Result ID	Reporting Name	LOINC®
FSCN4	SCN4A (Myotonia) DNA Sequencing	Not Provided