

Overview

Method Name

Immunofixation Electrophoresis

NY State Available

Yes

Specimen

Specimen Type

CSF

Specimen Required

Specimen Type: CSF

Container/Tube: sterile screw cap container

Specimen Volume: 3 mL

Collection Instructions: Collect 3 mL of spinal fluid (CSF) in a sterile screw cap container. Ship refrigerated.

Reject Due To

Hemolysis NA

Lipemia NA

Icterus NA

Other

Specimen Minimum Volume

1.5 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	7 days	
	Frozen	180 days	
	Ambient		

Clinical & Interpretive

Reference Values

No abnormal bands are present on immunofixation.

Performance**PDF Report**

Referral

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes**Test Classification****CPT Code Information**

86335

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FIXCF	Immunofixation, CSF	13174-8

Result ID	Test Result Name	Result LOINC Value
FIXCF	Immunofixation, CSF	13174-8