

**Overview****Method Name**

Blocking Radioimmunoassay (RIA) after acid-alcohol extraction

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain, red-top tube(s). Separate within 1 hour of collection, freeze immediately. Send 0.5 mL serum frozen.

Note: 1. Serum gel tube is okay, but must pour off into a plastic screw cap vial and freeze.  
2. Minimum volume does not permit for repeat analysis

**Specimen Minimum Volume**

0.1 mLNOTE: Minimum volume does not allow for repeat analysis.

**Reject Due To**

Hemolysis	NA
Lipemia	NA
Icteric	NA
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	200 days	

**Clinical & Interpretive****Reference Values**

---

<u>Age</u>	<u>Range (ng/mL)</u>	<u>Mean</u>
Prepubertal	258 – 882	570
Pubertal	273 – 892	583
Adults	333 – 967	650

## Performance

### PDF Report

No

### Day(s) Performed

Tuesday

### Report Available

3 to 14 days

### Performing Laboratory Location

Esoterix Endocrinology

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information

83519

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIGF2	Insulin-Like Growth Factor II	2485-1

Result ID	Test Result Name	Result LOINC® Value
Z2324	Insulin-Like Growth Factor II	2485-1