

Overview**Method Name**

Spectrophotometry

NY State Available

Yes

Specimen**Specimen Type**

Urine

Specimen Required

Container/Tube: Plastic urine container

Specimen Volume: 2 mL

Collection Instructions:

1. Collect random urine without preservative.
2. Ship 2 mL urine in a sterile screw capped plastic container ship ambient.

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	Preservatives

Specimen Minimum Volume

1 mL

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Ambient (preferred)	5 days	
	Frozen	21 days	

	Refrigerated	5 days	
--	--------------	--------	--

Clinical & Interpretive**Reference Values**

Adult: < or = 4 U/L

Performance**PDF Report**

No

Specimen Retention Time**Performing Laboratory Location**

Quest Diagnostics Nichols Institute

Fees & Codes**Test Classification****CPT Code Information**

83690

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FLIPR	Lipase, Urine	10888-6

Result ID	Reporting Name	LOINC®
Z1559	Lipase, Urine	10888-6