

Reporting Title: Nuts Allergen Profile

Performing Location Rochester

Ordering Guidance:

For a listing of allergens available for testing, see [Allergens - Immunoglobulin E \(IgE\) Antibodies](#) in Special Instructions

Specimen Requirements:

Container/Tube:

Preferred: Red top

Acceptable: Serum gel

Specimen Volume: 0.5 mL for every 5 allergens requested

Forms:

[If not ordering electronically, complete, print, and send an Allergen Test Request](#) (T236) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	90 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC
ALM	Almond, IgE	Numeric	kU/L	6019-4
CCNT	Coconut, IgE	Numeric	kU/L	6081-4
PEAN	Peanut, IgE	Numeric	kU/L	6206-7
PEC	Pecan-Food, IgE	Numeric	kU/L	6208-3
SESA	Sesame Seed, IgE	Numeric	kU/L	6242-2

Supplemental Report:

No

Components:

Test Id	Reporting Name	Available Separately	Always Performed
ALM	Almond, IgE	Yes	Yes
SESA	Sesame Seed, IgE	Yes	Yes
PEAN	Peanut, IgE	Yes	Yes
PEC	Pecan-Food, IgE	Yes	Yes
CCNT	Coconut, IgE	Yes	Yes

CPT Code Information:

86003 x 5

Reference Values:

For allergens ALM, SESA, PEC and CCNT:

Class	IgE kU/L	Interpretation
0	<0.35	Negative
1	0.35-0.69	Equivocal

2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive
6	> or =100	Strongly positive

Reference values apply to all ages.

For Allergen PEAN:

Class	IgE kU/L	Interpretation
0	<0.10	Negative
0/1	0.10-0.34	Borderline/Equivocal
1	0.35-0.69	Equivocal
2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive
6	> or =100	Strongly positive

Reference values apply to all ages.