

---

**Reporting Title:** Helminthosporium sativum IgE**Performing Location** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Supplemental Report:**

No

**CPT Code Information:**

86003

**Reference Values:**

&lt;0.35 kU/L