

Reporting Title: DNA Extract and Hold

Performing Location: Rochester

Specimen Requirements:

Specimen must arrive within 7 days of collection. Collect and package specimen as close to shipping time as possible.

Submit only 1 of the following specimens:

Specimen Type: Blood

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD)

Specimen Volume: 4 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send specimen in original tube. **Do not** aliquot.
3. Label specimen as blood.

Specimen Type: Bone marrow

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD)

Specimen Volume: 2 mL

Collection Instructions:

1. Invert several times to mix bone marrow.
2. Send specimen in original tube. **Do not** aliquot.
3. Label specimen as bone marrow.

Forms:

1. [Hematopathology Patient Information](#) (T676) in Special Instructions

2. If not ordering electronically, complete, print, and send a [Hematopathology/Cytogenetics Test Request](#) (T726) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	PURPLE OR PINK TOP/EDTA
	Refrigerated	7 days	PURPLE OR PINK TOP/EDTA

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
EXHD	MP034	Specimen Type	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC
MP034	Specimen Type	Alphanumeric		31208-2

41995	DNA Extract and Hold Result	Alphanumeric		No LOINC Needed
38305	DNA Extraction	Alphanumeric		No LOINC Needed

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:**Reference Values:**

Not applicable