

**Reporting Title:** Hazelnut Component Panel

**Performing Location:** Eurofins Viracor

**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC
Z5641	Hazelnut Component rCor a 1	Alphanumeric		69421-6
Z5642	Class	Alphanumeric		81995-3
Z5643	Hazelnut Component rCor a 8	Alphanumeric		58753-5
Z5644	Class	Alphanumeric		81993-8
Z5645	Hazelnut Component rCor a 9	Alphanumeric		65765-0
Z5646	Class	Alphanumeric		82002-7
Z5647	Hazelnut Component rCor a 14	Alphanumeric		81788-2
Z5648	Class	Alphanumeric		82569-5

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86008 x 4

**Reference Values:**

Class IgE	(kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal/Borderline
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	High Positive
4	17.50-49.99	Very High Positive

5	50.00-99.99	Very High Positive
6	>99.99	Very High Positive