
Reporting Title: Angiotensin Convert Enzyme CSF**Performing Location:** ARUP Laboratories**Specimen Requirements:****Container/Tube:** Sterile vial**Specimen Volume:** 1 mL**Collection Instructions:**

1. Transfer 1 mL of spinal fluid to a sterile vial.
2. Send frozen.

Note: Gadolinium contrast agents have been reported to inhibit ACE activity. Therefore, CSF containing gadolinium-based contrast agents should not be submitted to the laboratory for evaluation.

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	180 days	
	Refrigerated	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FACEC	Angiotensin Convert Enzyme CSF	Alphanumeric		12480-0

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82164

Reference Values:

0.0-2.5 U/L