

**Reporting Title:** Rhodotorula IgE  
**Performing Location:** Eurofins Viracor

**Specimen Requirements:**

**Specimen Type:** Serum

**Container/Tube:** Red or SST

**Specimen Volume:** 0.5 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z4036	Rhodotorula IgE	Alphanumeric		7665-3
Z4037	CLASS	Alphanumeric		102572-5

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86003

**Reference Values:**

<0.35 kU/L