

# **Test Definition: RPRT1**

Rapid Plasma Reagin (RPR) Screen with Reflex to Titer, Serum

## **Reporting Title:** RPR Screen w/ Reflex to Titer,S **Performing Location:** Rochester

## **Ordering Guidance:**

This test is for monitoring response to therapy in patients treated for syphilis infection. This test should not be used as a primary diagnostic approach for syphilis. To screen for undiagnosed syphilis infection, order SYPH1 / Syphilis IgG with Reflex, Enzyme Immunoassay, Serum.

### **Specimen Requirements:**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.5 mL Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

#### Forms:

If not ordering electronically, complete, print, and send <u>Infectious Disease Serology Test Request</u> (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
616863	RPR Screen w/ Reflex to Titer,S	Alphanumeric		20507-0

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

# **CPT Code Information:**

86592 86593-Rapid Plasma Reagin Titer (if appropriate)

# Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
RPRT2	RPR Titer,S	1	86593	No	No



# **Test Definition: RPRT1**

Rapid Plasma Reagin (RPR) Screen with Reflex to Titer, Serum

## **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
RPRT2	616865	RPR Titer,S	Alphanumeric		31147-2

## **Reference Values:**

Negative

Reference values apply to all ages.