

Reporting Title: Hepatosplenomegaly Panel, P

Performing Location: Rochester

Ordering Guidance:

This test **should not be used** for monitoring of patients with confirmed diagnoses. If testing requested is for monitoring purposes, see:

CTXP / Cerebrotendinous Xanthomatosis, Plasma

GPSYP / Glucopsychosine, Plasma

OXNP / Oxysterols, Plasma

Specimen Requirements:

- Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)
- Collection Container/Tube:**
- Preferred:** Lavender top (EDTA)
- Acceptable:** Green top (sodium heparin, lithium heparin), yellow top (ACD B)
- Submission Container/Tube:** Plastic vial
- Specimen Volume:** 0.3 mL
- Collection Instructions:**
1. Centrifuge at 4 degrees C, if possible
 2. Aliquot plasma into plastic vial. **Do not disturb or transfer the buffy coat layer.**
 3. Send frozen

Forms:

[If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request](#) (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma	Frozen	65 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
601542	Interpretation (HSMP)	Alphanumeric		59462-2
601536	Cholestane-3beta,5alpha,6beta-triol	Numeric	nmol/mL	92755-8
601537	7-Ketocholesterol	Numeric	nmol/mL	92764-0
601538	Lyso-sphingomyelin	Numeric	nmol/mL	92747-5
601539	Glucopsychosine	Numeric	nmol/mL	92750-9
601540	7a-hydroxy-4-cholesten-3-one	Numeric	nmol/mL	92761-6
601541	7a,12a-dihydroxycholest-4-en-3-one	Numeric	nmol/mL	92758-2
601543	Reviewed By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:
82542

Reference Values:

Cholestane -3-beta, 5-alpha, 6-beta-troil
Cutoff: < or =0.070 nmol/mL

7-Ketocholesterol
Cutoff: < or =0.100 nmol/mL

Lyso-sphingomyelin Cutoff: < or =0.100 nmol/mL

Glucopsychosine Cutoff: < or =0.003 nmol/mL

7-Alpha-hydroxy-4-cholesten-3-one Cutoff: < or =0.300 nmol/mL

7-Alpha,12-alpha-dihydroxycholest-4-en-3-one
Cutoff: < or =0.100 nmol/mL