



Test Definition: FBILM

Biliary Tract Malignancy-Cytology, FISH, Varies

Reporting Title: Biliary Tract Malignancy-Cyto/FISH

Performing Location: Mayo Clinic Laboratories - Rochester Main Campus

Specimen Requirements:

Supplies:

PreservCyt Vial (T536)

CytoLyt Solution (T564)

Specimen Type: Bile duct brushing, bile duct aspirate, hepatobiliary brushing, or hepatobiliary aspirate

Container/Tube: Separate ThinPrep vial containing 20 mL PreservCyt or CytoLyt solution for each specimen

Specimen Volume: Entire collection

Collection Instructions: Label with site specimen was collected from (eg, right hepatic duct or common bile duct).

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|------|-------------------|
| Varies | Refrigerated (preferred) | | |
| | Ambient | | |

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description | Type | Reportable |
|---------|-------------|----------------------|------------|------------|
| FBILM | CY070 | Collection Procedure | Plain Text | Yes |
| FBILM | CY042 | Source | Plain Text | Yes |
| FBILM | CY043 | Clinical History | Plain Text | Yes |
| FBILM | CY044 | Fixative | Plain Text | Yes |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|------------------------------------|--------------|------|-----------------|
| CY070 | Collection Procedure | Alphanumeric | | 33724-6 |
| CY042 | Source | Alphanumeric | | 22633-2 |
| CY043 | Clinical History | Alphanumeric | | 22636-5 |
| CY044 | Fixative | Alphanumeric | | 8100-0 |
| 71816 | Case Number | Alphanumeric | | 80398-1 |
| 71291 | Interpretation | Alphanumeric | | 69965-2 |
| 71292 | Participated in the Interpretation | Alphanumeric | | No LOINC Needed |
| 71293 | Report electronically signed by | Alphanumeric | | 19139-5 |
| 71294 | Addendum | Alphanumeric | | 35265-8 |
| 71295 | Gross Description | Alphanumeric | | 22634-0 |
| 71570 | Disclaimer | Alphanumeric | | 62364-5 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88112
88377 (if appropriate)

Reflex Tests:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Available Separately |
|---------|--------------------------------|-----------|----------|------------------|----------------------|
| BILMB | Biliary Tract Malignancy, FISH | 1 | 88377 | No | No |
| BILMC | Biliary Tract Malignancy, FISH | | | No | No |
| BILMD | Biliary Tract Malignancy, FISH | 1 | 88377 | No | No |
| BILME | Biliary Tract Malignancy, FISH | 1 | 88377 | No | No |
| BILMF | Biliary Tract Malignancy, FISH | 1 | 88377 | No | No |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name | Type | Unit | LOINC® |
|---------|-----------|---------------------|--------------|------|---------|
| BILMA | 71864 | Result Summary | Alphanumeric | | 50397-9 |
| BILMA | 71865 | Result | Alphanumeric | | 95231-7 |
| BILMA | 71866 | Interpretation | Alphanumeric | | 69965-2 |
| BILMA | 71867 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILMA | 71868 | Specimen | Alphanumeric | | 31208-2 |
| BILMA | 71869 | Source | Alphanumeric | | 31208-2 |
| BILMA | 71870 | Method | Alphanumeric | | 85069-3 |
| BILMA | 71871 | Released By | Alphanumeric | | 18771-6 |
| BILMB | 71873 | Result Summary | Alphanumeric | | 50397-9 |
| BILMB | 71874 | Result | Alphanumeric | | 95231-7 |
| BILMB | 71875 | Interpretation | Alphanumeric | | 69965-2 |
| BILMB | 71876 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILMB | 71877 | Specimen | Alphanumeric | | 31208-2 |
| BILMB | 71878 | Source | Alphanumeric | | 31208-2 |
| BILMB | 71879 | Method | Alphanumeric | | 85069-3 |
| BILMB | 71880 | Released By | Alphanumeric | | 18771-6 |
| BILMC | 71882 | Result Summary | Alphanumeric | | 50397-9 |
| BILMC | 71883 | Result | Alphanumeric | | 95231-7 |
| BILMC | 71884 | Interpretation | Alphanumeric | | 69965-2 |
| BILMC | 71885 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILMC | 71886 | Specimen | Alphanumeric | | 31208-2 |
| BILMC | 71887 | Source | Alphanumeric | | 31208-2 |
| BILMC | 71888 | Method | Alphanumeric | | 85069-3 |
| BILMC | 71889 | Released By | Alphanumeric | | 18771-6 |
| BILMD | 71891 | Result Summary | Alphanumeric | | 50397-9 |
| BILMD | 71892 | Result | Alphanumeric | | 95231-7 |

| | | | | | |
|-------|-------|---------------------|--------------|--|---------|
| BILMD | 71893 | Interpretation | Alphanumeric | | 69965-2 |
| BILMD | 71894 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILMD | 71895 | Specimen | Alphanumeric | | 31208-2 |
| BILMD | 71896 | Source | Alphanumeric | | 31208-2 |
| BILMD | 71897 | Method | Alphanumeric | | 85069-3 |
| BILMD | 71898 | Released By | Alphanumeric | | 18771-6 |
| BILMF | 71909 | Result Summary | Alphanumeric | | 50397-9 |
| BILMF | 71910 | Result | Alphanumeric | | 95231-7 |
| BILMF | 71912 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILMF | 71911 | Interpretation | Alphanumeric | | 69965-2 |
| BILMF | 71913 | Specimen | Alphanumeric | | 31208-2 |
| BILMF | 71914 | Source | Alphanumeric | | 31208-2 |
| BILMF | 71915 | Method | Alphanumeric | | 85069-3 |
| BILMF | 71916 | Released By | Alphanumeric | | 18771-6 |
| BILME | 71900 | Result Summary | Alphanumeric | | 50397-9 |
| BILME | 71901 | Result | Alphanumeric | | 95231-7 |
| BILME | 71902 | Interpretation | Alphanumeric | | 69965-2 |
| BILME | 71903 | Reason for Referral | Alphanumeric | | 42349-1 |
| BILME | 71904 | Specimen | Alphanumeric | | 31208-2 |
| BILME | 71905 | Source | Alphanumeric | | 31208-2 |
| BILME | 71906 | Method | Alphanumeric | | 85069-3 |
| BILME | 71907 | Released By | Alphanumeric | | 18771-6 |

Reference Values:

Negative for malignancy.