

Reporting Title: IGG Subtypes Immunofluorescence
Performing Location: Mayo Clinical Laboratories- Rochester Main Campus

Ordering Guidance:
If additional interpretation/analysis is needed, request PATHC / Pathology Consultation along with this test and send the corresponding renal pathology light microscopy and immunofluorescence (IF) slides (or IF images on a CD), electron microscopy images (prints or CD), and the pathology report.

Shipping Instructions:
1. Advise shipping frozen specimens (unstained slides or tissue block) in Styrofoam transportation coolers filled with dry ice to ensure specimens are received at required specimen stability temperature.
2. Attach the green "Attention Pathology" address label (T498) to the outside of the transport container before putting into the courier mailer.

Necessary Information:
A pathology/diagnostic report is required.

Specimen Requirements:
Submit only 1 of the following specimens:

Preferred:
Specimen Type: Unstained slides (unfixed)
Source: Kidney tissue
Slides: 4 Slides
Collection Instructions: 4 frozen tissue unstained positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick, **centered on the slide**, and submitted on dry ice.

Acceptable:
Specimen Type: Unfixed tissue block (O.C.T)
Source: Kidney tissue
Specimen Volume: Entire specimen
Collection Instructions:
1. Embed in O.C.T. compound.
2. Freeze specimen and ship on dry ice.

Acceptable:
Specimen Type: Wet tissue
Source: Kidney tissue
Supplies: Renal Biopsy Kit (T231)
Container/Tube: Renal Biopsy Kit, Zeus/Michel's
Specimen Volume: Entire specimen
Collection Instructions:
1. Collect specimens according to the instructions in [Renal Biopsy Preparation Instructions](#).
2. If standard immunoglobulin and complement immunofluorescence has already been performed, submit the residual frozen tissue (must contain glomeruli) on dry ice.

Forms:

1. [Renal Biopsy Patient Information](#)
2. If not ordering electronically, complete, print, and send a [Renal Diagnostics Test Request](#) (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Special	Frozen		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
BA0271	Interpretation	Alphanumeric		50595-8
BA0272	Participated in the Interpretation	Alphanumeric		No LOINC Needed
BA0276	Addendum	Alphanumeric		35265-8
BA0275	Report electronically signed by	Alphanumeric		19139-5
71618	Disclaimer	Alphanumeric		62364-5
71850	Case Number	Alphanumeric		80398-1
BA0274	Gross Description	Alphanumeric		22634-0
BA0273	Material Received	Alphanumeric		22633-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88346-Primary IF
88350-If additional IF

Reference Values:

An interpretive report will be provided.