

Reporting Title: NuSwab Vaginitis (VG)
Performing Location: LabCorp Burlington

Specimen Requirements:

Collection Container/Tube: APTIMA Multitest, vaginal, or unisex swab

Specimen Volume: One swab

Collection Instructions: Collect vaginal fluid sample using the Gen-Probe Aptima swab by contacting the swab to the lower third of the vaginal wall and rotating the swab for 10 to 30 seconds to absorb fluid. Immediately place the swab into the transport tube and carefully break the swab shaft against the side of the tube. Tightly screw on the cap. Submit one vaginal swab in APTIMA Multitest, vaginal, or unisex swab. Ship refrigerate.

Specimen Type	Temperature	Time	Special Container
Swab	Ambient	30 days	
	Refrigerated (preferred)	30 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z4735	Atopobium vaginae	Alphanumeric		69565-0
Z4736	BVAB 2	Alphanumeric		69566-8
Z4737	Megasphaera 1	Alphanumeric		69567-6
Z4738	Candida albicans, NAA	Alphanumeric		69562-7
Z4739	Candida glabrata, NAA	Alphanumeric		69563-5
Z4740	Trich vag by NAA	Alphanumeric		62461-9

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87801
87798 x 3
87661

Reference Values:

Candida albicans, NAA: Negative
Candida glabrata, NAA: Negative
Trich vag by NAA: Negative