



Test Definition: FARP

Autoimmune Retinopathy Panel by
Immunoblot (ARP)

Reporting Title: Autoimmune Retinopathy Panel (ARP)

Performing Location: Ocular Immunology Laboratory OHSU

Specimen Requirements:

Submit only one of the following specimens:

Serum:

Draw blood in a plain, red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 5 mL serum refrigerated in a plastic vial.

Plasma:

Draw blood in a lavender-top (EDTA) tube(s). Spin down and send 5 mL EDTA plasma refrigerated in a plastic vial.

Complete and submit with specimen:

1. Completed OHSU Ocular request form
2. Clinical history
3. Referring physician information (name & phone number)

-NOTE: Without this information, testing cannot be completed.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FARP	Autoimmune Retinopathy Panel (ARP)	Alphanumeric		Not Provided

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

Referral

CPT Code Information:

84182 x 8

Reference Values:

A final report will be provided.