



Test Definition: FEBGP

Exotic Bird Panel IgG

Reporting Title: Exotic Bird Panel IgG

Performing Location: Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5631	Canary Feathers IgG	Alphanumeric		58916-8
Z5632	Finch Feathers IgG	Alphanumeric		Not Provided
Z5633	Parrot Aust Drp	Alphanumeric		Not Provided
Z5634	Pigeon Feathers IgG	Alphanumeric		58912-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86001 x4

Reference Values:

Canary Feathers IgG	<17.0
Finch Feathers IgG	<12.0
Parrot Australian (Budgerigar) Droppings IgG	<14.0
Pigeon Feathers IgG	<22.0