



# Test Definition: HAPT

Haptoglobin, Serum

**Reporting Title:** Haptoglobin, S

**Performing Location:** Mayo Clinic Jacksonville Clinical Lab

**Specimen Requirements:**

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL Serum

**Collection Instructions:** Centrifuge and aliquot serum into a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	14 days	
	Frozen	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
HAPT	Haptoglobin, S	Numeric	mg/dL	46127-7

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83010

**Reference Values:**

30-200 mg/dL