

Test Definition: LAGU

Legionella Antigen, Random, Urine

Reporting Title: Legionella Ag, U **Performing Location:** Jacksonville

Ordering Guidance:

This assay has been verified using urine specimens only. For serum specimens, order SLEG / Legionella pneumophila (Legionnaires Disease), Antibody, Serum.

Other specimen types (eg, plasma or body fluids) that may contain Legionella antigen have not been verified for testing.

Specimen Requirements:

Supplies: Urine Tubes, 10 mL (T068) **Container/Tube:** Plastic, 10-mL urine tube

Specimen Volume: 0.5 mL **Collection Instructions:**

- 1. Collect a random urine specimen.
- 2. No preservative.
- 3. Excessively bloody or very turbid specimens containing protein, cells, or particulates will be canceled. They can inhibit the function of the test.
- 4. Centrifuging to remove particulates is not approved.
- 5. Specimens with any dyes or unnatural color are not acceptable and will be canceled.

Forms:

If not ordering electronically, complete, print, and send 1 of the following with the specimen:

- -Infectious Disease Serology Test Request (T916)
- -Microbiology Test Request (T244)

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|----------|-------------------|
| Urine | Ambient | 24 hours | |
| | Refrigerated (preferred) | 7 days | |
| | Frozen | 14 days | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|------------------|--------------|------|---------|
| 81268 | Legionella Ag, U | Alphanumeric | | 32781-7 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87899



Test Definition: LAGU

Legionella Antigen, Random, Urine

| _ | - | | | | |
|---|-------------|-------|-----|------------|-------|
| D | ∧ f∧ | ror | 200 | V 2 | lues: |
| Г | ele | :: =: | ıce | va | iues. |

Negative