



Test Definition: FLUNP

Influenza Virus Type A and Type B, and Respiratory Syncytial Virus (RSV), Molecular Detection, PCR, Nasopharyngeal Swab

Reporting Title: Influenza A/B and RSV, PCR, NP Swab

Performing Location: Mayo Clinic Jacksonville Clinical Lab

Specimen Requirements:

Specimen Type: Nasopharyngeal swab

Container/Tube: Sterile container with viral transport media

Collection Instructions:

1. Collect specimen by swabbing back and forth over mucosa surface to maximize recovery of cells.
2. Swab must be placed in viral transport media (for example, M4-RT, M4, or M4 media).

Forms:

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|--------|-------------------|
| Swab | Refrigerated (preferred) | 7 days | |
| | Frozen | 7 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|----------------------------------|--------------|------|---------|
| SS019 | Specimen Source | Alphanumeric | | 31208-2 |
| 35975 | Influenza A, PCR | Alphanumeric | | 76078-5 |
| 35976 | Influenza B, PCR | Alphanumeric | | 76080-1 |
| 35977 | Respiratory Syncytial Virus, PCR | Alphanumeric | | 76089-2 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87631

Reference Values:

Negative