

## **Test Definition: MELAI**

Melaleuca leucadendron, IgE, Serum

### Reporting Title: Melaleuca Leucadendron, IgE

Performing Location: Mayo Clinic Laboratories - Rochester Superior Drive

#### Ordering Guidance:

For a listing of allergens available for testing, see <u>Allergens - Immunoglobulin E (IgE) Antibodies</u>.

#### **Specimen Requirements:**

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.5 mL for every 5 allergens requested Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

#### Forms:

If not ordering electronically, complete, print, and send an <u>Allergen Test Request</u> (T236) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	90 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
MELAI	Melaleuca Leucadendron, IgE	Numeric	kU/L	6171-3

LOINC® and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

# **CPT Code Information:** 86003

#### **Reference Values:**

Class	lgE kU/L	Interpretation
0	<0.10	Negative
0/1	0.10-0.34	Borderline/equivocal
1	0.35-0.69	Equivocal
2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive



## **Test Definition: MELAI**

Melaleuca leucadendron, IgE, Serum

6	> or =100	Strongly positive
U U	× 01 ±00	

Reference values apply to all ages.

Concentrations of 0.70 kU/L or more (class 2 and above) will flag as abnormally high.