

Reporting Title: Diazepam and Nordiazepam, S

Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube: Red top (Serum gel/SST are not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions: Centrifuge and aliquot serum into plastic vial.

Forms:

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum Red     | Ambient                  | 14 days |                   |
|               | Refrigerated (preferred) | 14 days |                   |
|               | Frozen                   | 14 days |                   |

Result Codes:

| Result ID | Reporting Name           | Type         | Unit  | LOINC®  |
|-----------|--------------------------|--------------|-------|---------|
| 8629      | Diazepam                 | Numeric      | ng/mL | 3548-5  |
| 2475      | Nordiazepam              | Numeric      | ng/mL | 3537-8  |
| 2459      | Diazepam and Nordiazepam | Alphanumeric | ng/mL | 16757-7 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80299

Reference Values:

Therapeutic concentrations

Diazepam and Nordiazepam: 200-2,500 ng/mL