

**Reporting Title:** Glucose, Random, S  
**Performing Location:** Mayo Clinic Laboratories - Rochester Main Campus

**Specimen Requirements:**

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.5 mL

**Collection Instructions:**

- 1. Serum gel tubes should be centrifuged within 2 hours of collection.
- 2. Red-top tubes should be centrifuged, and the serum aliquoted into a plastic vial within 2 hours of collection.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	30 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
GLURA	Glucose, Random, S	Numeric	mg/dL	2345-7

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

82947

**Reference Values:**

- 0-11 months: not established
- > or =1 year: 70-140 mg/dL