



Test Definition: FFURO

Furosemide (Lasix)

Reporting Title: Furosemide (Lasix)

Performing Location: Medtox Laboratories, Inc.

Specimen Requirements:

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 1 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 1 mL of serum refrigerated in a plastic vial.

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|----------|-------------------|
| Varies | Refrigerated (preferred) | 7 days | |
| | Ambient | 72 hours | |
| | Frozen | 180 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|----------------|--------------|------|--------|
| Z1189 | Furosemide | Alphanumeric | | 3659-0 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80299

Reference Values:

Expected serum furosemide concentration in patients on usual daily dosages: Up to 5.0 ug/mL

Toxic: greater than 50.0 ug/mL