

Reporting Title: Calcitonin, S
Performing Location: Rochester

Specimen Requirements:

Patient Preparation: For 12 hours before specimen collection, patient should not take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

- 1. After collection, immediately place specimen on ice.
- 2. Refrigerate specimen during centrifugation and immediately transfer serum to a plastic vial.

Forms:

If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	8 hours	
	Refrigerated	24 hours	
	Frozen (preferred)	90 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
CATN	Calcitonin, S	Numeric	pg/mL	1992-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82308

Reference Values:

Pediatric

1 month: < or =34 pg/mL

2 months: < or =31 pg/mL

3 months: < or =28 pg/mL

4 months: < or =26 pg/mL

5 months: < or =24 pg/mL

6 months: < or =22 pg/mL
7 months: < or =20 pg/mL
8 months: < or =19.0 pg/mL
9 months: < or =17.0 pg/mL
10 months: < or =16.0 pg/mL
11 months: < or =15.0 pg/mL
12-14 months: < or =14.0 pg/mL
15-17 months: < or =12.0 pg/mL
18-20 months: < or =10.0 pg/mL
21-23 months: < or =9.0 pg/mL
2 years: < or =8.0 pg/mL
3-9 years: < or =7.0 pg/mL
10-15 years: < or =6.0 pg/mL
16 years: < or =5.0 pg/mL

Adults

17 years and older:

Males: < or =14.3 pg/mL

Females: < or =7.6 pg/mL

For International System of Units (SI) conversion for Reference Values, see
www.mayocliniclabs.com/order-tests/si-unit-conversion.html.