

Reporting Title: Prostate-Specific Ag Diagnostic, S**Performing Location:** Mayo Clinic Laboratories - Rochester Superior Drive**Necessary Information:**

Include patient's age.

Specimen Requirements:**Supplies:** Sarstedt Aliquot Tube, 5mL (T914)**Collection Container/Tube:****Preferred:** Serum gel**Acceptable:** Red top**Submission Container/Tube:** Plastic vial**Specimen Volume:** 0.6 mL serum**Collection Instructions:**

1. For serum gel tubes, centrifuge within 2 hours of collection.
2. For red top tubes, centrifuge and aliquot serum into a plastic vial within 2 hours of collection.

Additional Information: Free prostate-specific antigen (PSA) can only be added on within 72 hours of performing total PSA. Specimen must have been shipped frozen.**Forms:**If not ordering electronically, complete, print, and send an [Oncology Test Request](#) (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	180 days	
	Refrigerated	5 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
PSA	Prostate-Specific Ag Diagnostic, S	Numeric	ng/mL	83112-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84153

Reference Values:

Males:

<40 years: < or =2.0 ng/mL

40-49 years: < or =2.5 ng/mL

50-59 years: < or =3.5 ng/mL
60-69 years: < or =4.5 ng/mL
70-79 years: < or =6.5 ng/mL
> or =80 years: < or =7.2 ng/mL

Females: Not applicable