



Complete and return this form with your specimen. Specimens received without this information may be rejected.

Instructions for the Patient

- **Do not** brush teeth, eat, drink, or take any oral medication for at least **60 minutes** before collecting specimen.
- **Rinse** mouth thoroughly with water for 1–5 seconds, then wait **10 minutes** before collecting specimen.
- **Collect specimen at the time instructed by your healthcare professional.** Collection immediately prior to going to bed is most common, but collection immediately after waking is appropriate in some circumstances.
- **Document** Patient ID or Birth Date on the container label of the swab tube.
- **Record** the exact collection **date and time** on **BOTH** the container label (swab tube) and in the box below:

Collection Date (mm-dd-yyyy)

Collection Time (hh:mm)

☐ am
☐ pm

Barcode
Label



Cap


Swab

Insert

Tube



1. Remove cap to expose swab. **Do not** remove the plastic insert from the tube. **Do not** use fingers to touch the swab.



2. Place swab directly into the mouth by tipping the tube. **Do not** let fingers touch the swab. Keep swab in mouth for **2 minutes**. Roll swab in mouth. **Do not** chew.



3. Spit swab back into tube. **Do not** let fingers touch the swab.



4. Replace the cap. Make sure it is placed on tightly.



5. Record the **exact** collection date and time (include **am** or **pm**) on barcode label **and** on this form.

Return the tube (with swab) and this form to your healthcare provider within **24 to 48 hours**.