

Overview**Method Name**

Enzymatic/Spectrophotometric

NY State Available

Yes

Specimen**Specimen Type**

Fecal

Specimen Required

Collect 1 gm random stool in sterile leak proof container, ship refrigerate.

Note:

Dietary restrictions: Patients receiving pancreatic enzymes should discontinue taking the enzymes at least 5 days before the collection of the stool sample.

Specimen Minimum Volume

0.5 gram

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Fecal	Refrigerated (preferred)	14 days	
	Frozen	30 days	

Clinical and Interpretive**Reference Values**

2.3 - 51.4 U/g

Performance

PDF Report

No

Day(s) and Time(s) Test Performed

Wednesday

Analytic Time

2 days

Maximum Laboratory Time

4 - 18 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FCHYS	Chymotrypsin, Stool	25375-7

Result ID	Test Result Name	Result LOINC Value
FCHYS	Chymotrypsin, Stool	25375-7