

**Overview**
**Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen**
**Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Hemolysis:	Mild OK; Gross reject
Thawing:	Warm OK; Cold OK
Lipemia:	Mild OK; Gross reject
Icterus:	Mild OK; Gross reject
Other:	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Clinical and Interpretive**
**Reference Values**

&lt;0.35 kU/L

**Interpretation**

Class IgE (kU/L) Comment

0 &lt;0.10 Negative

0/1 0.10  $\leq$  0.34 Equivocal/Borderline

1 0.35  $\leq$  0.69 Low Positive

2 0.70  $\leq$  3.49 Moderate Positive

3 3.50  $\leq$  17.49 High Positive

4 17.50  $\leq$  49.99 Very High Positive

5 50.00  $\leq$  99.99 Very High Positive

6 >99.99 Very High Positive

## Performance

### PDF Report

No

### Day(s) and Time(s) Test Performed

Monday through Friday

### Analytic Time

1 - 2 days

### Maximum Laboratory Time

3 - 4 days

### Performing Laboratory Location

Viracor Eurofins

## Fees and Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

86003

### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FIGE	Fig IgE	11173-2



Result ID	Test Result Name	Result LOINC Value
Z3870	Fig IgE	11173-2
Z3871	CLASS	21270-4