

Overview**Method Name**

ImmunoCAP FEIA

NY State Available

Yes

Specimen**Specimen Type**

Serum

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Clinical and Interpretive**Reference Values**

<0.35 kU/L

Interpretation

Class IgE (kU/L) Comment

0 <0.10 Negative

0/1 0.10 - 0.34 Equivocal/Borderline

1 0.35 - 0.69 Low Positive

2 0.70 - 3.49 Moderate Positive

3 3.50 - 17.49 High Positive

4 17.50 - 49.99 Very High Positive

5 50.00 - 99.99 Very High Positive

6 >99.99 Very High Positive

Performance

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Friday

Analytic Time

1 - 2 days

Maximum Laboratory Time

3 - 6 days

Performing Laboratory Location

Eurofins Viracor

Fees and Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86003

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FOAKE	Oak Live (Quercus virginiana) IgE	6164-8

Result ID	Test Result Name	Result LOINC Value
Z4118	Oak Live (<i>Quercus virginiana</i>) IgE	6164-8
Z4119	CLASS	15880-8