

Overview**Method Name**

Hemolysis

NY State Available

Yes

Specimen**Specimen Type**

Serum Red

Specimen Required

Draw blood in a plain red-top tube. Maintain specimen at 37 degrees C until serum is separated from cells. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

1.5 mL

Reject Due To

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	Separator or gel tube

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

Clinical and Interpretive**Reference Values**

Negative

Performance**PDF Report**

No

Day(s) and Time(s) Test Performed

Monday through Friday

Analytic Time

1 - 3 days

Maximum Laboratory Time

3 - 7 days

Performing Laboratory Location

ARUP Laboratories

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

CPT Code Information

86940, 86941

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FDLS	Donath Landsteiner	18287-3

Result ID	Test Result Name	Result LOINC Value
FDLS	Donath Landsteiner	18287-3